


ROCKER (NB-104)



Features

- No vibration, No malfunction
- Safe using for Petri-dish, culture bottle and flask
- Microprocessor control, Using 10 ~ 35 rpm, setting up to 3hour



Specification



Items	Unit	NB-104
Speed range	rpm	10 to 35 rpm
Angle range	°	±7° from horizontal
Time range	hour	Up to 3 hours or continuous
Platform size	mm	290(W) x 200(D)mm
Capacity		96 well micro-plate x 4 ea
Dimensions	mm	290(W) x 210(D) x 145(H)mm
Power	V/Hz	220V / 60Hz
Weight	kg	5kg